

Sensoronics

Model SP0-500 Wrist Oximeter



Technical Parameters

The one to be worn comfortably on the patients wrist for easy spot-check of blood oxygen saturation and pulse rate.

Can be used for monitoring daily activities and overnight sleep studies,

The compact design is intended to give the patient greater freedom while improving patient compliance.

Compact and light weight design, fits comfortably on the patient wrist

Simply plug in the sensor

For spot checking short and long term monitoring applications

High accuracy and repetition

Friendly operation menu for the function setting

Sensoronics

Incorporated

5114 Pt. Fosdick Dr. NW, Suite E70
Gig Harbor, WA 98355

Phone 877-677-0816

Fax 877-677-0136

www.sensoronics.com